



MicroTaker/BX<sup>®</sup> for Buffer  
Exchange Applications


# Overview

We need a way to automate the isolating and washing of nanoparticles or biomolecules post harvesting.  
Also known as 'Buffer Exchange' Applications.


# Workflow Overview

Nanoparticle Isolating and Washing Workflow		
Step	Action	Est. Time
M-1	Load 2 racks of 24 1 dram vials with fabricated nanoparticles onto the robot deck along with on filtration plate and initiate the program	N/A
A-1	The system uncaps all the vessels and places the caps on a separate rack	552 Sec.
A-2	The UpLyft rack tips angling the vials to that the most sample can be aspirated via pipette and transferred to the associated position on the filtration plate	16 Sec.
A-3	This is repeated for the first 24 nanoparticle samples	368 Sec.
A-4	DMSO system solvent is selected ant the line washed and primed	25 Sec.
A-5	X mL of DMSO is added to each well of the filter plate in the shaker	50 Sec.
A-6	Shaker shakes vials for 5 minutes and then a gradual vacuum is applied removing the DMSO int the waste line	05 Sec.
A-7	Steps A-5 - A-6 is repeated two more times	165 Sec.
A-8	Steps A-4 - A-7 is repeated with a high volatile solvent.	245 Sec.
M-2	Operator is notified that the first rack is complete and that a new filter rack and cap rack need to be added and the program resumed.	N/A
A-9	Process is completed for the next 24 nanoparticles	1426 Sec.
<b>Total Est. Time</b>		<b>2852 Sec.</b>
<b>Total Est. Time</b>		<b>47.5 Min.</b>

# Recommended Configuration



# Vacuum Filtration with Shaking



# Recommended Platform



- Standard and pipette tip cannula
- Capping/uncapping station for 1-dram vials
- Pick and place for 1-dram vials
- Custom filtration-based washing/solvent exchange station