







1. 3D Discovery Tools

Our most popular product line helps researchers and project managers from cancer discovery to High-Throughput-Screening.





Prof. Taniguchi and Dr. Takebe, Department of Regenerative Medicine, Yokohama City University.

2. Spheroid Generators


The device can be placed on 6-well plate. Hundreds of spheroids can be generated by just seeding cell suspension in the device. Harvesting the spheroids is so easy.

Hundreds of spheroids can be generated by just seeding cell suspension in the device.









D=350µm (Product number: MPc 350 6)
D=500µm (Product number: MPc 500 6)
D=700µm (Product number: MPc 700 6)

Place the device into medium in another dish



Product list

1. 3D Discovery Tools

SEM Micrograph	Surface	Plate format	Product number
	Non-adherent surface for spheroid culture	6-well	SQ 200 100 NA 6
		24-well	SQ 200 100 NA 24
		96-well	SQ 200 100 NA 96
		384-well	SQ 200 100 NA 384
	Adherent surface	6-well	SQ 200 100 6
		24-well	SQ 200 100 24
		96-well	SQ 200 100 96
		384-well	SQ 200 100 384
	Non-adherent surface for spheroid culture	6-well	RB 500 400 NA 6
		24-well	RB 500 400 NA 24
		96-well	RB 500 400 NA 96
	Adherent surface	6-well	RB 500 400 6
		24-well	RB 500 400 24
		96-well	RB 500 400 96
	Non-adherent surface for spheroid culture	6-well	RB 900 700 NA 6
		24-well	RB 900 700 NA 24
	Adherent surface	6-well	RB 900 700 6
		24-well	RB 900 700 24
	Adherent surface	24-well	SQ 200 50 24
		96-well	SQ 200 50 96
		384-well	SQ 200 50 384
	Adherent surface	24-well	SQ 100 50 24
		96-well	SQ 100 50 96
		384-well	SQ 100 50 384
	Adherent surface	24-well	CG 85 25 20 24
		96-well	CG 85 25 20 96
		384-well	CG 85 25 20 384

2. Spheroid Generators



Set	Product number
6 devices/plate	MPc 350 6
6 devices/plate	MPc 500 6
6 devices/plate	MPc 700 6

Price and availability in your country, please contact with us.

Elplasia.3d@kuraray.com , <http://www.elplasia.com>