

Sample collection and processing labels

Well and microplate labels

Brady's well and microplate labels can withstand varying temperature processes and have a thin profile, short height and long width to fit the necessary information on the plate.

- B-488 – White polyester with moderate chemical resistance and short-term exposure to xylene, DMSO and ethanol
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers




Figure 2




Figure 25




Figure 82




Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges								
Fig. 2	M-81-488	B-488	Black on white	1.900 (48.3)	0.250 (6.4)	1	160	In cartridge
BradyPrinter M611 label rolls								
Fig. 82	M61-2-7425	B-7425	White	0.250 (6.4)	2.000 (50.8)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls								
Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	2	1,500	B30-R6000
Fig. 2	B33-151-490	B-490	White	2.000 (50.8)	0.250 (6.4)	1	2,500	B30-R4300
Fig. 2	B33-151-499	B-499	White	2.000 (50.8)	0.250 (6.4)	1	1,500	B30-R6000
Fig. 2	B33-151-7425	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	B30-R4300
Fig. 2	B33-53-7425	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	B30-R4300
BradyJet J2000 Printer label rolls								
Fig. 2	J20-151-7425J	B-7425J	White	2.000 (50.8)	0.250 (6.4)	1	2,000	J20-CMY
Fig. 2	J20-53-7425J	B-7425J	White	2.000 (50.8)	0.500 (12.7)	1	2,000	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls								
Fig. 2	THF-151-499-1.5-SC	B-499	White	2.000 (50.8)	0.250 (6.4)	1	1,500	R4913
Fig. 2	THF-151-7425-2-SC	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R6413
Fig. 2	THF-53-7425-2-SC	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R6413
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls								
Fig. 2	THF-15-488-2.5	B-488	White	2.000 (50.8)	0.250 (6.4)	1	2,500	R4302 / IP-R4302
Fig. 2	THF-151-499-3	B-499	White	2.000 (50.8)	0.250 (6.4)	1	3,000	R6000 / IP-R6000
Fig. 2	THF-151-7425-2	B-7425	White	2.000 (50.8)	0.250 (6.4)	1	2,000	R4302 / IP-R4302
Fig. 2	THF-53-7425-2	B-7425	White	2.000 (50.8)	0.500 (12.7)	1	2,000	R4302 / IP-R4302

Centrifuge and conical tube labels

Designed with ultra-thin material that prevents vials from jamming in centrifuges or racks. These labels withstand harsh temperatures and chemicals without falling off, fading or smearing.

- B-461 – All-white, ultra-thin polyester label. Also available in a self-laminating format
- B-427 – Durable all-white vinyl label is good for wrapping and conformability. Must be wrapped upon itself when used in liquid nitrogen
- B-432 – Clear polyester withstands exposure to solvents and exhibits good adhesion to metal and other surfaces
- B-490 – White polyester ideal for flat frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 0.125 inch
- B-492 – Ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature tubes and vials. Must have small overlap if applied to already frozen containers
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges										
15 mL	Fig. 33b	M21-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges										
15 mL	Fig. 83	M-95-492	B-492	Black on white	0.600 (15.2)	1.625 (41.3)	0.375 (9.5) dia.	1	200	In cartridge
	Fig. 2	M-156-492	B-492	Black on white	1.000 (25.4)	0.900 (22.9)	-	1	180	In cartridge
	Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 82	M-78-432-CL-BK	B-432	Black on clear	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 82	M-78-492	B-492	Black on white	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 82	M-128-499	B-499	Black on white	1.000 (25.4)	1.900 (48.3)	-	1	100	In cartridge
	Fig. 33	MC-1000-499	B-499	Black on white	1.000 (25.4)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-1000-461	B-461	Black on white	1.000 (25.4)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-1250-461	B-461	Black on white	1.250 (31.8)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
50 mL	Fig. 82	M-113-490	B-490	Black on white	1.500 (38.1)	2.200 (55.9)	-	1	140	In cartridge
	Fig. 82	M-114-490	B-490	Black on white	1.500 (38.1)	3.750 (95.3)	-	1	80	In cartridge

Centrifuge and conical tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Centrifuge and conical tube labels (continued)




Figure 2




Figure 3




Figure 5a




Figure 33




Figure 82

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BradyPrinter M611 label rolls										
15 mL	Fig. 82	M61-156-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 82	M61-203-490	B-490	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R4310
	Fig. 82	M61-203-492	B-492	White	1.000 (25.4)	1.750 (44.5)	-	1	250	M61-R6410
	Fig. 82	M61-204-492	B-492	White	1.000 (25.4)	1.900 (48.3)	-	1	250	M61-R6410
	Fig. 33	M61C-1000-499	B-499	White	1.000 (25.4)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-1000-461	B-461	All White	1.000 (25.4)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 33	M61C-1250-461	B-461	All White	1.250 (31.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
50 mL	Fig. 2	PTL-78-490	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	250	M61-R4310
	Fig. 82	PTL-113-490	B-490	White	1.500 (38.1)	2.200 (55.9)	-	1	100	M61-R4310
	Fig. 82	PTL-114-490	B-490	White	1.500 (38.1)	3.875 (98.4)	-	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
15 mL	Fig. 2	B33-179-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	1,900	B30-R6400
	Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-136-499	B-499	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	B30-R6000
	Fig. 2	B33-136-7425	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	B30-R4300
	Fig. 2	B33-258-7425	B-7425	White	1.500 (38.1)	1.750 (44.5)	-	1	900	B30-R4300
	Fig. 2	B33-155-490	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	B30-R4300
	Fig. 2	B33-156-490	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-156-492	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	B30-R6400
50 mL	Fig. 2	B33-137-499	B-499	White	2.000 (50.8)	1.000 (25.4)	-	1	1,100	B30-R6000
	Fig. 2	B33-137-7425	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-157-490	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,000	B30-R4300
	Fig. 2	B33-158-490	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	1,000	B30-R4300
BradyJet J2000 Printer label rolls										
15 mL	Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY
	Fig. 3	J20C-1125-7425J	B-7425J	White	1.125 (28.6)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
	Fig. 2	J20-136-7425J	B-7425J	White	1.500 (38.1)	0.750 (19.1)	-	1	1,500	J20-CMY
	Fig. 2	J20-258-7425J	B-7425J	White	1.500 (38.1)	1.750 (44.5)	-	1	900	J20-CMY
	Fig. 3	J20C-2250-7425J	B-7425J	White	2.250 (57.2)	100 ft. (30.5 m)	-	1	1 roll (cont.)	J20-CMY
50 mL	Fig. 2	J20-137-7425J	B-7425J	White	2.000 (50.8)	1.000 (25.4)	-	1	1,500	J20-CMY
	Fig. 2	J20-18-7425J	B-7425J	White	3.000 (76.2)	1.000 (25.4)	-	1	1,000	J20-CMY

Centrifuge and conical tube labels continue on next page.




Sample collection and processing labels

Centrifuge and conical tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
15 mL	Fig. 2	THT-179-492-2.5-SC	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	2,500	R6413
	Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
	Fig. 2	THT-136-7425-2-SC	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	2,000	R6413
	Fig. 2	THT-155-490-1.5-SC	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	1,500	R4313
	Fig. 2	THT-258-7425-2-SC	B-7425	White	1.750 (44.5)	1.500 (38.1)	-	1	2,000	R6413
	Fig. 2	THT-156-490-1.5-SC	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R4313
	Fig. 2	THT-156-492-1.5-SC	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	1,500	R6413
15 mL top	Fig. 5a	THT-143-492-1.5-SC	B-492	White	0.375 (9.5) dia.	-	-	4	1,500	R6413
50 mL	Fig. 2	THT-137-7425-2-SC	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R6413
	Fig. 2	THT-157-490-1.5-SC	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	1,500	R4313
	Fig. 2	THT-158-490-1.5-SC	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
15 mL	Fig. 2	THT-135-427-3-AW	B-427	All White	0.750 (19.1)	1.500 (38.1)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-179-492-3	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-136-7425-2	B-7425	White	1.500 (38.1)	0.750 (19.1)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-155-490-3	B-490	White	1.625 (41.3)	0.600 (15.2)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-258-7425-2	B-7425	White	1.750 (44.5)	1.500 (38.1)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-156-490-3	B-490	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-156-492-3	B-492	White	1.900 (48.3)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
50 mL	Fig. 2	THT-137-7425-2	B-7425	White	2.000 (50.8)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-137-427-3-AW	B-427	All White	2.000 (50.8)	1.000 (25.4)	-	1	3,000	R4302 / IP-R4302
	Fig. 2	THT-157-490-3	B-490	White	2.200 (55.9)	1.500 (38.1)	-	1	3,000	R4300 / IP-R4300
	Fig. 2	THT-158-490-3	B-490	White	3.750 (95.3)	1.500 (38.1)	-	1	3,000	R4307 / IP-R4307

Self-laminating centrifuge and conical tube labels

- B-461 – Ultra-thin polyester with clear overlaminate tail that protects print and allows liquid levels to be visible



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges										
15 mL	Fig. 15	M-129-461	B-461	Black on white	2.250 (57.2)	1.500 (38.1)	1.250 (31.7)	1	140	In cartridge
	Fig. 15	M-119-461	B-461	Black on white	3.750 (95.3)	1.500 (38.1)	1.250 (31.7)	1	80	In cartridge
BradyPrinter M611 label rolls										
15 mL	Fig. 14	PTL-77-461	B-461	White	3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	1	250	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
15 mL	Fig. 14c	B33-276-461	B-461	White	2.625 (66.7)	1.000 (25.4)	1.000 (25.4)	2	1,500	B30-R4300
	Fig. 15	B33-134-461	B-461	White	3.500 (88.9)	1.500 (38.1)	1.250 (31.8)	1	1,000	B30-R4300
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
15 mL	Fig. 14c	THT-139-461-2	B-461	White	1.000 (25.4)	5.250 (133.4)	1.500 (38.1)	2	2,000	R4300 / IP-R4300
	Fig. 15	THT-134-461-3	B-461	White	1.500 (38.1)	3.500 (88.9)	1.250 (31.8)	1	3,000	R4307 / IP-R4307


*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels

Designed to fit the small sizes of microcentrifuge and PCR tubes, these labels withstand harsh temperatures and chemicals without falling off, fading or smearing. Once in the microcentrifuge, vial top labels allow you to quickly extract specific vials.

- B-430 – Clear high-performance polyester label
- B-461 – All-white, ultra-thin polyester label. Also available in a self-laminating format
- B-490 – White polyester is ideal for flat frozen surfaces. Adheres to frozen tubes and vials when wrapped around itself by at least 0.125 inch
- B-492 – Ultra-thin white polyester with superior chemical and abrasion resistance for labeling frosted, frozen or room-temperature tubes and vials. Must have small overlap if applied to already frozen containers
- B-494 – Glossy polyester with a printed color bar that allows for quick visual differentiation between samples. Withstands liquid nitrogen, a variety of chemicals and more. Label dimensions listed below include the color bar area. Vial tops are solid color
- B-499 – White nylon cloth with permanent adhesive for easy handling and conformability
- B-7425 / B-7425J – Economical white polypropylene label with moderate chemical resistance that is ideal for labeling dry, room-temperature specimen containers. Must be applied to containers prior to being frozen



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®21-LAB Label Printer cartridges										
2 mL	Fig. 33b	M21-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-430	B-430	Black on clear	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-7425	B-7425	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-750-461	B-461	Black on white	0.750 (19.1)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
0.5 mL	Fig. 33b	M21-250-430	B-430	Black on clear	0.250 (6.4)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-499	B-499	Black on white	0.375 (9.5)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-430	B-430	Black on clear	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-7425	B-7425	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-375-461	B-461	Black on white	0.375 (9.5)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-430	B-430	Black on clear	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-7425	B-7425	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-461	B-461	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33b	M21-500-461	B-461	Black on white	0.500 (12.7)	21 ft. (6.4 m)	-	1	1 roll (cont.)	In cartridge
BMP®51 Label Printer cartridges										
0.5 mL	Fig. 82	M-130-499	B-499	Black on white	0.375 (9.5)	0.825 (21.0)	-	1	220	In cartridge
	Fig. 33	MC-375-499	B-499	Black on white	0.375 (9.5)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-375-461	B-461	Black on white	0.375 (9.5)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-499	B-499	Black on white	0.500 (12.7)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-492	B-492	Black on white	0.500 (12.7)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 33	MC-500-461	B-461	Black on white	0.500 (12.7)	25 ft. (7.6 m)	-	1	1 roll (cont.)	In cartridge

Microcentrifuge and PCR tube labels continue on next page.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.5 mL	Fig. 82	M-131-492	B-492	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-BL	B-494	Black on blue / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-GN	B-494	Black on green / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-OR	B-494	Black on orange / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-PK	B-494	Black on pink / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-RD	B-494	Black on red / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-494-YL	B-494	Black on yellow / white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-131-499	B-499	Black on white	0.500 (12.7)	1.000 (25.4)	-	1	180	In cartridge
	Fig. 82	M-124-492	B-492	Black on white	0.500 (12.7)	1.650 (41.9)	-	1	180	In cartridge
	Fig. 24	M-118-492	B-492	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-GN	B-494	Black on green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-YL	B-494	Black on red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-118-494-RD	B-494	Black on yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	250	In cartridge
	Fig. 24	M-118-499	B-499	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	240	In cartridge
	Fig. 24	M-118-7425	B-7425	Black on white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-BL	B-494	Black on blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-GN	B-494	Black on green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-OR	B-494	Black on orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-PK	B-494	Black on pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-RD	B-494	Black on red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-494-YL	B-494	Black on yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-120-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	240	In cartridge
	Fig. 24	M-117-492	B-492	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 24	M-117-499	B-499	Black on white	1.000 (25.4)	0.500 (12.7)	0.438 (11.1) dia.	2	240	In cartridge
	Fig. 2	M-112-490	B-490	Black on white	1.250 (31.8)	0.375 (9.5)	-	1	360	In cartridge
	Fig. 2	M-112-492	B-492	Black on white	1.250 (31.8)	0.375 (9.5)	-	1	360	In cartridge
2 mL	Fig. 33	MC-750-499	B-499	Black on white	0.750 (19.1)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 2	M-156-492	B-492	Black on white	1.000 (25.4)	0.900 (22.9)	-	1	180	In cartridge
	Fig. 2	M-49-494-BL	B-494	Black on blue / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-GN	B-494	Black on green / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-OR	B-494	Black on orange / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-PK	B-494	Black on pink / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-RD	B-494	Black on red / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
	Fig. 2	M-49-494-YL	B-494	Black on yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	260	In cartridge
Any top	Fig. 33	MC-1000-499	B-499	Black on white	1.000 (25.4)	16 ft. (4.8 m)	-	1	1 roll (cont.)	In cartridge
	Fig. 8	M-82-499	B-499	Black on white	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-492	B-492	Black on white	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-OR-BK	B-499	Black on orange	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-RD-BK	B-499	Black on red	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-82-499-YL-BK	B-499	Black on yellow	0.375 (9.5) dia.	-	-	1	240	In cartridge
	Fig. 8	M-83-499	B-499	Black on white	0.500 (12.7) dia.	-	-	1	240	In cartridge
	Fig. 8	M-83-492	B-492	Black on white	0.500 (12.7) dia.	-	-	1	240	In cartridge
	Fig. 8	M-83-499-OR-BK	B-499	Black on orange	0.500 (12.7) dia.	-	-	1	240	In cartridge
	Fig. 8	M-83-499-RD-BK	B-499	Black on red	0.500 (12.7) dia.	-	-	1	240	In cartridge
0.5 mL	Fig. 8	M-83-499-YL-BK	B-499	Black on yellow	0.500 (12.7) dia.	-	-	1	240	In cartridge
	BradyPrinter M611 label rolls									
	Fig. 33	M61C-375-499	B-499	White	0.375 (9.5)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-375-461	B-461	White	0.375 (9.5)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
	Fig. 33	M61C-500-499	B-499	White	0.500 (12.7)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 33	M61C-500-492	B-492	White	0.500 (12.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R6410
	Fig. 33	M61C-500-461	B-461	White	0.500 (12.7)	50 ft. (15.2 m)	-	1	1 roll (cont.)	M61-R4310
0.5 mL	Fig. 2	M61-130-499	B-499	White	0.825 (21.0)	0.375 (9.5)	-	1	500	M61-R6210

Microcentrifuge and PCR tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)




Figure 1




Figure 2




Figure 5a




Figure 8




Figure 24




Figure 25




Figure 33




Figure 82

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.5 mL	Fig. 82	M61-47-499	B-499	White	0.500 (12.7)	1.000 (25.4)	-	1	500	M61-R6210
	Fig. 82	M61-124-492	B-492	White	0.500 (12.7)	1.650 (41.9)	-	1	500	M61-R6410
	Fig. 24	M61-117-492	B-492	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6410
	Fig. 24	PTL-117-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	500	M61-R6210
	Fig. 24	M61-98-492	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	M61-98-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6410
	Fig. 24	PTL-98-499	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R6210
	Fig. 24	M61-98-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	500	M61-R4310
	Fig. 2	M61-17-492	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	250	M61-R6410
	Fig. 2	M61-17-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R6410
	Fig. 2	M61-17-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	500	M61-R4310
	Fig. 24	M61-120-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6410
	Fig. 24	M61-120-499	B-499	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	250	M61-R6210
	Fig. 2	M61-199-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	250	M61-R4310
	Fig. 2	PTL-112-490	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	500	M61-R4310
2 mL	Fig. 33	M61C-750-499	B-499	White	0.750 (19.1)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 82	M61-156-492	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-18-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	1	250	M61-R4310
	Fig. 2	M61-19-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 2	M61-19-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	250	M61-R6410
	Fig. 33	M61C-1000-499	B-499	White	1.000 (25.4)	30 ft. (9.1 m)	-	1	1 roll (cont.)	M61-R6210
	Fig. 8	M61-82-492	B-492	White	0.375 (9.5) dia.	-	-	1	500	M61-R6410
	Fig. 8	PTL-82-499	B-499	White	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	PTL-82-499-YL	B-499	Yellow	0.375 (9.5) dia.	-	-	1	500	M61-R6210
Any top	Fig. 8	M61-82-499-OR	B-499	Orange	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-82-499-RD	B-499	Red	0.375 (9.5) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-492	B-492	White	0.500 (12.7) dia.	-	-	1	500	M61-R6410
	Fig. 8	PTL-83-499	B-499	White	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	PTL-83-499-YL	B-499	Yellow	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-499-OR	B-499	Orange	0.500 (12.7) dia.	-	-	1	500	M61-R6210
	Fig. 8	M61-83-499-RD	B-499	Red	0.500 (12.7) dia.	-	-	1	500	M61-R6210
BradyPrinter i3300 and BBP®33 Printer label rolls										
0.5 mL	Fig. 1	B33-68-490	B-490	White	0.500 (12.7)	0.750 (19.1)	-	4	2,500	B30-R4300
	Fig. 1	B33-68-499	B-499	White	0.500 (12.7)	0.750 (19.1)	-	4	1,500	B30-R6000
	Fig. 2	B33-257-7425	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	B30-R4300
	Fig. 2	B33-125-490	B-490	White	0.900 (22.9)	0.500 (12.7)	-	1	1,650	B30-R4300
	Fig. 2	B33-125-499	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	B30-R6000
	Fig. 24	B33-126-490	B-490	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	2,500	B30-R4300
	Fig. 24	B33-126-499	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	2	1,500	B30-R4300

Microcentrifuge and PCR tube labels continue on next page.

*Tube sizes may vary by brand.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
1.5 mL	Fig. 25	B33-152-492	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	B30-R6400
	Fig. 25	B33-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6400
	Fig. 25	B33-152-499	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	1,500	B30-R6000
	Fig. 25	B33-152-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	B30-R4300
	Fig. 24	B33-163-492	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 2	B33-163-499	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R6000
	Fig. 24	B33-249-7425	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	B30-R4300
	Fig. 2	B33-59-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	B30-R4300
	Fig. 24	B33-181-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
	Fig. 24	B33-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	B30-R6400
2 mL	Fig. 24	B33-182-492	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R6400
	Fig. 24	B33-251-7425	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	B30-R4300
	Fig. 2	B33-154-490	B-490	White	1.250 (31.8)	0.375 (9.5)	-	1	1,500	B30-R4300
	Fig. 2	B33-256-7425	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	B30-R4300
	Fig. 25	B33-88-7425	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	B30-R4300
	Fig. 2	B33-179-492	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	1,900	B30-R6400
	Fig. 2	B33-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
Any top	Fig. 2	B33-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R6400
	Fig. 2	B33-179-7425	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	B30-R4300
	Fig. 5a	B33-143-492	B-492	White	0.375 (9.5)	-	-	4	1,500	B30-R6400
	Fig. 5a	B33-143-499	B-499	White	0.375 (9.5)	-	-	4	1,500	B30-R6000
	Fig. 5a	B33-144-490	B-490	White	0.500 (12.7) dia.	-	-	3	2,500	B30-R4300
BradyJet J2000 Printer label rolls	Fig. 5a	B33-144-492	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6400
	Fig. 5a	B33-144-499	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	B30-R6000
	Fig. 2	J20-257-7425J	B-7425J	White	0.825 (21.0)	0.250 (6.4)	-	1	1,500	J20-CMY
	Fig. 25	J20-152-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	J20-CMY
	Fig. 24	J20-249-7425J	B-7425J	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	2	1,500	J20-CMY
	Fig. 2	J20-59-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	J20-CMY
	Fig. 24	J20-251-7425J	B-7425J	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	2	1,500	J20-CMY
	Fig. 2	J20-256-7425J	B-7425J	White	1.300 (33.0)	0.600 (15.2)	-	1	1,500	J20-CMY
2 mL	Fig. 25	J20-88-7425J	B-7425J	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	J20-CMY
	Fig. 2	J20-179-7425J	B-7425J	White	1.000 (25.4)	1.000 (25.4)	-	1	1,500	J20-CMY
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
0.5 mL	Fig. 2	THT-257-7425-2-SC	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R6413
	Fig. 2	THT-125-499-1.5-SC	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	1,500	R4913
1.5 mL	Fig. 25	THT-152-492-2.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	2,500	R6413
	Fig. 25	THT-152-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	1,500	R4913
	Fig. 25	THT-152-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R6413
	Fig. 24	THT-249-492-1.5-SC	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-249-499-1.5-SC	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	1,500	R4913
	Fig. 24	THT-249-7425-2-SC	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R6413
	Fig. 2	THT-59-492-2.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2,500	R6413
	Fig. 2	THT-59-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R6413
	Fig. 24	THT-250-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	1,500	R6413
	Fig. 24	THT-251-492-1.5-SC	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	1,500	R6413
	Fig. 24	THT-251-7425-2-SC	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R6413
	Fig. 2	THT-256-7425-2-SC	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R6413

*Tube sizes may vary by brand.

Microcentrifuge and PCR tube labels continue on next page.

Sample collection and processing labels

Microcentrifuge and PCR tube labels (continued)




Figure 2




Figure 5a




Figure 24




Figure 25

Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Secondary label C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
2 mL	Fig. 25	THT-88-7425-2-SC	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R6413
	Fig. 2	THT-179-492-2.5-SC	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	2,500	R6413
	Fig. 2	THT-179-7425-2-SC	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R6413
Any top	Fig. 5a	THT-143-492-1.5-SC	B-492	White	0.375 (9.5) dia.	-	-	4	1,500	R6413
	Fig. 5a	THT-144-492-1.5-SC	B-492	White	0.500 (12.7) dia.	-	-	3	1,500	R6413
	Fig. 5a	THT-144-499-1.5-SC	B-499	White	0.500 (12.7) dia.	-	-	3	1,500	R4913
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
0.5 mL	Fig. 2	THT-257-7425-2	B-7425	White	0.825 (21.0)	0.250 (6.4)	-	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-125-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	-	1	3,000	R6000 / IP-R6000
	Fig. 24	THT-126-499-3	B-499	White	0.900 (22.9)	0.500 (12.7)	0.500 (12.7) dia.	1	3,000	R6000 / IP-R6000
1.5 mL	Fig. 25	THT-152-492-3	B-492	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-BL	B-494	Blue / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-GN	B-494	Green / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-OR	B-494	Orange / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-PK	B-494	Pink / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-RD	B-494	Red / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-494-YL	B-494	Yellow / white	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6402 / IP-R6402
	Fig. 25	THT-152-499-3	B-499	White	1.000 (25.4)	0.375 (9.5)	-	2	3,000	R6000 / IP-R6000
	Fig. 25	THT-152-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	-	2	2,000	R4302 / IP-R4302
	Fig. 24	THT-163-492-3	B-492	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-163-499-3	B-499	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	3,000	R6000 / IP-R6000
	Fig. 24	THT-249-7425-2	B-7425	White	1.000 (25.4)	0.375 (9.5)	0.375 (9.5) dia.	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-59-492-10	B-492	White	1.000 (25.4)	0.500 (12.7)	-	1	2 rolls of 5,000	R6400 / IP-R6400
	Fig. 2	THT-59-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	-	1	2,000	R4302 / IP-R4302
	Fig. 24	THT-181-492-3	B-492	White	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6400 / IP-R6400
	Fig. 24	THT-250-494-BL	B-494	Blue / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-GN	B-494	Green / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-OR	B-494	Orange / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-PK	B-494	Pink / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-RD	B-494	Red / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-250-494-YL	B-494	Yellow / white	1.000 (25.4)	0.500 (12.7)	0.375 (9.5) dia.	1	3,000	R6402 / IP-R6402
	Fig. 24	THT-182-492-3	B-492	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	3,000	R6400 / IP-R6400
	Fig. 24	THT-251-7425-2	B-7425	White	1.000 (25.4)	0.500 (12.7)	0.440 (11.2) dia.	1	2,000	R4302 / IP-R4302
	Fig. 2	THT-256-7425-2	B-7425	White	1.300 (33.0)	0.600 (15.2)	-	1	2,000	R4302 / IP-R4302
2 mL	Fig. 2	THT-180-492-3	B-492	White	0.900 (22.9)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 25	THT-88-7425-2	B-7425	White	1.000 (25.4)	0.750 (19.1)	-	2	2,000	R4302 / IP-R4302
	Fig. 2	THT-179-492-3	B-492	White	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6400 / IP-R6400
	Fig. 2	THT-179-494-BL	B-494	Blue / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-GN	B-494	Green / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-OR	B-494	Orange / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-PK	B-494	Pink / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-RD	B-494	Red / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-494-YL	B-494	Yellow / white	1.000 (25.4)	1.000 (25.4)	-	1	3,000	R6402 / IP-R6402
	Fig. 2	THT-179-7425-2	B-7425	White	1.000 (25.4)	1.000 (25.4)	-	1	2,000	R4302 / IP-R4302
Any top	Fig. 5a	THT-143-492-3	B-492	White	0.375 (9.5) dia.	-	-	3	3,000	R6402 / IP-R6402
	Fig. 5a	THT-143-499-3	B-499	White	0.375 (9.5) dia.	-	-	4	3,000	R6000 / IP-R6000
	Fig. 5a	THT-144-492-3	B-492	White	0.500 (12.7) dia.	-	-	3	3,000	R6402 / IP-R6402
	Fig. 5a	THT-144-499-3	B-499	White	0.500 (12.7) dia.	-	-	3	3,000	R6000 / IP-R6000

Self-laminating microcentrifuge and PCR tube labels

- B-461 – Ultra-thin polyester with clear overlamine tail that protects print and allows liquid levels to be visible




Figure 14




Figure 14b




Figure 15



Fits tube size*	Figure	Catalog #	Material	Color	Label width A inch (mm)	Label length B inch (mm)	Print-on area C inch (mm)	Labels per row	Labels per roll	Recommended ribbon
BMP®51 Label Printer cartridges										
0.5 mL	Fig. 15	M-123-461	B-461	Black on white	1.125 (28.6)	0.375 (9.5)	0.750 (19.0)	1	360	In cartridge
	Fig. 15	M-122-461	B-461	Black on white	1.125 (28.6)	0.500 (12.7)	0.750 (19.0)	1	360	In cartridge
	Fig. 15	M-124-461	B-461	Black on white	1.650 (41.9)	0.500 (12.7)	0.750 (19.0)	1	180	In cartridge
	Fig. 15	M-72-461	B-461	Black on white	1.750 (44.5)	0.500 (12.7)	0.750 (19.0)	1	180	In cartridge
1.5 mL	Fig. 15	M-126-461	B-461	Black on white	1.800 (45.7)	0.600 (15.2)	1.000 (25.4)	1	180	In cartridge
	Fig. 15	M-86-461	B-461	Black on white	2.625 (66.7)	0.600 (15.2)	1.000 (25.4)	1	120	In cartridge
2 mL	Fig. 15	M-74-461	B-461	Black on white	2.200 (55.9)	0.500 (12.7)	0.750 (19.0)	1	140	In cartridge
BradyPrinter M611 label rolls										
0.5 mL	Fig. 14	M61-202-461	B-461	White	0.500 (12.7)	1.650 (41.9)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	M61-123-461	B-461	White	1.125 (28.6)	0.375 (9.5)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	M61-122-461	B-461	White	1.125 (28.6)	0.500 (12.7)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	PTL-76-461	B-461	White	1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	1	250	M61-R4310
	Fig. 15	PTL-73-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	500	M61-R4310
	Fig. 15	PTL-72-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	500	M61-R4310
1.5 mL	Fig. 14	M61-201-461	B-461	White	0.600 (15.2)	1.800 (45.7)	1.000 (25.4)	1	100	M61-R4310
	Fig. 15	PTL-108-461	B-461	White	1.500 (38.1)	0.187 (4.75)	0.700 (17.8)	1	750	M61-R4310
2 mL	Fig. 14	M61-203-461	B-461	White	1.000 (25.4)	1.750 (44.5)	1.000 (25.4)	1	100	M61-R4310
	Fig. 14	PTL-74-461	B-461	White	2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	1	100	M61-R4310
	Fig. 14	PTL-79-461	B-461	White	3.750 (95.3)	0.750 (19.1)	1.000 (25.4)	1	100	M61-R4310
BradyPrinter i3300 and BBP®33 Printer label rolls										
0.5 mL	Fig. 15	B33-123-461	B-461	White	1.800 (45.7)	0.375 (9.5)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 15	B33-122-461	B-461	White	1.800 (45.7)	0.500 (12.7)	0.750 (19.1)	1	1,500	B30-R4300
	Fig. 14b	B33-127-461	B-461	White	2.200 (55.9)	0.500 (12.7)	0.750 (19.1)	3	1,500	B30-R4300
1.5 mL	Fig. 14b	B33-133-461	B-461	White	1.625 (41.3)	0.600 (15.2)	1.000 (25.4)	3	1,500	B30-R4300
BBP®12 Printer and other industry printers using 1-inch small core label rolls										
0.5 mL	Fig. 15	THT-122-461-1.5-SC	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-127-461-1.5-SC	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	1,500	R4313
	Fig. 15	THT-123-461-1.5-SC	B-461	White	0.750 (19.1)	1.800 (45.7)	0.750 (19.1)	1	1,500	R4313
1.5 mL	Fig. 15	THT-133-461-1.5-SC	B-461	White	0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	1	1,500	R4313
	Fig. 15	THT-153-461-1.5-SC	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	1,500	R4313
BradyPrinter i5100, i7100, Brady IP Printers and other industry printers using 3-inch core label rolls										
0.5 mL	Fig. 15	THT-123-461-3	B-461	White	0.375 (9.5)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-122-461-3	B-461	White	0.500 (12.7)	1.800 (45.7)	0.750 (19.1)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-127-461-3	B-461	White	0.500 (12.7)	2.200 (55.9)	0.750 (19.1)	1	3,000	R4300 / IP-R4300
1.5 mL	Fig. 15	THT-133-461-3	B-461	White	0.600 (15.2)	1.625 (41.3)	1.000 (25.4)	1	3,000	R4302 / IP-R4302
	Fig. 15	THT-153-461-3	B-461	White	0.600 (15.2)	2.625 (66.7)	1.000 (25.4)	1	3,000	R4300 / IP-R4300

*Tube sizes may vary by brand.