




microPro 300


Features:


- Multiple Pipetting Modes
- Intuitive Touchscreen User Interface
- RingLock Tip Loading System
- Automatic Pipetting Depth Setting
- Automatic Tip Ejection
- Save Commonly Used Protocols
- 10 Pipetting Speeds for All Liquid Types
- 384 Well Microplate Pipetting



Specifications:

| | |
|------------|---|
| Volume | 5-300 μ L |
| Footprint | 19.5 cm x 21 cm (7.5" x 8.25") |
| Max Height | 31.5 cm (12.4") |
| Weight | 4 kg (10 lbs) |
| Precision | 15-30 μ L <2% C.V. 30-300 μ L <1% C.V. |





000.0 μ L

Reservoir


| | | | |
|-------------|---------------|-----------------|----------|
| Aspirate | 50.0 μ L | Reservoir | 5 |
| | Volume | Container | Speed |
| | 20.0 μ L | Air Gap | 5 |
| | 50.0 μ L | Rounded Bottom2 | 5 |
| | 120.0 μ L | none | 5 |
| + Add Mix | | | |
| New Program | | | |
| Purge | Save | Help | Settings |

Pipetting Modes:

- Simple Pipette with Mix
- Repeat Dispense
- Variable Dispense
- Dilute with Mix
- Multi Aspirate with Mix
- Custom Programs


96 and 384 Well Applications Include:

- Cell Culture: Media Exchange
- ELISA Assays
- Magnetic Bead Separation Assays
- 96 to 384 Daughter Plate Creation



We are a Massachusetts based company specializing in precision liquid handling products. Our mission is to create innovative solutions that simplify the jobs of researchers and lab technicians.

Quality, simplicity, and precision are the foundation for everything we make.



10 Technology Dr, Suite 40
Hudson, MA 01749
www.avidien.com