

(800) 344-1321 www.obtura.com

Obtura Spartan Ultrasonic Tips

Power Settings:

In general, start with a lower setting and only incrementally increase the power if necessary to safely accomplish the procedure.

*The power settings used are determined by the clinician's experience with ultrasonics and the power needed to safely complete the specific procedure. *Indicated power settings are intended for use with Obtura Spartan ultrasonic units. Setting may not be appropriate for other units.

*See Instruction Manual for Use.








gold plating that improves cutting



Allows enhanced vision for better access vs. high speed handpiece



BUC-2

Disk-like tip allows you to horizontally smooth pulp chamber floor

BUC-3

20.2mm

Extremely active tip for cutting apically into calcified canals and troughing around posts

BUC-2A

Smaller diameter allows for precise smoothing around corners

BUC-3A

Designed for troughing around broken instruments in the coronal half of the root

24.2mm

CPR-1

For vibrating posts and cores, and aiding in crown and bridge removal



Smaller taper sub-orifice tip used to remove obstructions in the canal



Longer than the CPR 6, but with the same flexibility, strength and cutting ability

CPR-2D

efficiency

CPR[®] Ultrasonic Tips

17.0mm

RETREATMENT

This versatile and durable tip provides enhanced energy for removal of core materials

17.0mm

CPR-5D

Smaller taper sub-orifice tip used to remove obstructions in the canal





17.0mm

CPR-6


CPR-3D

Smaller taper sub-orifice

Combines the flexibility and strength of a titanium alloy tip with the proven cutting ability of a CPR tip

Obturation **Ultrasonics** Surgery Irrigation Diagnostic:

The Gold Standard Ultrasonic Tip Guide



Ultrasonic Accessories

Tip Taxi

Designed to hold all of Obtura Spartan's endodontic ultrasonic instruments. Organizes tips and wrench, eliminating needless loss. Autoclavable.



Endo File Adapter

The Endo File Adapter is designed to function on most brands of Piezo-Electric type dental ultrasonic power units with standard M3x.6 threading. The adapter can accommodate K Type file sizes 15-40 once the handle or latch hub has been removed.

MTS Ultrasonic Unit

The MTS is the premier multi-task system for all procedures such as endodontics, periodontics, and restorative dentistry. When using magnification and Obtura Spartan Endodontic's specialized instruments, this complete system allows you to perform to your exact standards.



Tip Wrench

Specifically designed for all Obtura Spartan tips including: BUC Access Refinement, CPR, SL Micro-restorative, KiS Microsurgical, File Adapter and Scaling Tips.

BUC Ultrasonic Tips

930-041	BUC 1
930-048	BUC 1A
930-042	BUC 2
930-046	BUC 2A
930-043	BUC 3
930-049	BUC 3A
930-050	Set of 6 BUC Tips

CPR Ultrasonic Tips

30-011	CPR 1D
30-012	CPR 2D
30-013	CPR 3D
30-014	CPR 4D
30-015	CPR 5D
30-016	CPR 6 Titanium
30-017	CPR 7 Titanium
30-018	CPR 8 Titanium
30-115	CPR 1-5 Tip Set
30-608	CPR 6-8 Tip Set

KiS Ultrasonic Tips

931-011	KiS 1D
931-012	KiS 2D
931-013	KiS 3D
931-014	KiS 4D
931-015	KiS 5D
931-016	KiS 6D
931-116	KiS 1-6 Tip Se
952-611	KiS 1D2
952-613	KiS 3D2
952-614	KiS 1D4
952-615	KiS 5D2

PST Ultrasonic Tips

941-101	PST 1
941-102	PST 2
941-103	PST 3
941-104	PST 4
941-105	PST 5
941-106	PST 6
941-107	PST 7
941-108	PST 8
941-109	Set of 8 PST Tips

SL Ultrasonic Tips

0-031	SL 1A
0-032	SL 2A
0-033	SL 3A
0-034	SL 4A
0-135	Set of 4 SL Tips

Spartan MTS Ultrasonic Unit

42-100 Piezo-Electric MTS 115V43-710 Autoclavable Endo Handpiece

Ultrasonic Tip Accessories

-300 Tip Adapter Male/Male		
-001 Wrench Scaler Tip	-001	Wrench Scaler Tip

