

sigvaris

Lower extremity formulary guide

**A guide to selecting
compression solutions**



To Order:
1.800.322.7744, ext. 1



Online Dealer Portal:
sigvaris-online.com








sigvaris.com

Chronic venous disease product options

Medical compression is a firmly established treatment for venous and lymphatic conditions, as well as for venous symptoms in healthy people. Robust evidence is published regarding the use of medical compression stockings for management of venous symptoms, prevention and treatment of venous edema and management of venous ulcers (*Rabe et al., 2018*).*

◆ Day Wear ◆ Night Wear


	
Style Sheer	Style Soft Opaque
Fashionable, sheer medical graduated compression hosiery	Graduated compression with lasting microfiber softness, advanced moisture management, and opaque colors

Available Compression Strengths (measured in mmHg)					15-20	20-30	30-40	15-20	20-30	30-40
Chronic Venous Insufficiency	CEAP Classification**		Diagnosis	Aim of Compression						
	C0		No visible or palpable signs of venous disease but high risk from pregnancy or hereditary disease	Prevention of swelling, heavy legs, tenderness, pain, and DVT	◆			◆		
	C1		Telangiectasia or reticular veins (spider veins)	Prevention of swelling, heavy legs, tenderness and pain	◆	◆		◆	◆	
	C2		Varicose Veins	Vein support, prevention of swelling, heavy legs, tenderness, and pain	◆	◆	◆	◆	◆	◆
	C3		Edema	Edema reduction		◆	◆		◆	◆
	C4		3 subclasses: a) pigmentation changes (hemosiderin staining) or eczema b) Lipodermatosclerosis or atrophie blanche c) Corona Phlebectatica	Prevention of ulceration		◆	◆		◆	◆
	C5		Healed venous ulcer	Prevention of ulcer recurrence		◆	◆		◆	◆
	C6		Active venous ulcer	Ulcer closure, pain relief, vein support						
Other	Phlebo-Lymphedema		Accumulation of fluid, protein and other elements in the tissue/ interstitial spaces	Reduction of edema						

* Rabe et al. Indications for medical compression stockings in venous and lymphatic disorders: An evidence-based consensus statement. Phlebology, 2018, Vol. 33(3) 163-184

** CEAP Classification System and Reporting Standard Revision, 2020

															
Essential Cotton		Essential Opaque		Style Microfiber & Style Microfiber Patterns			Secure			Compreflex®			Coolflex™		
Breathable and durable supima cotton stockings and socks designed for everyday wear and sensitive skin		From petite to plus, this garment offers the largest assortment of sizes and styles of graduated compression stockings		Soft and durable sock for everyday wear, neutral and exciting colors and patterns for men and women			Graduated compression and containment hybrid garment, most containing circular knit garments			Gradient, short-stretch compression; each garment offers Accutab® compression measuring devices			Graduated compression with visual indicator of compression level; breathable mesh for cool and comfortable fit		
20-30	30-40	20-30	30-40	15-20	20-30	30-40	20-30	30-40	40-50	20-30	30-40	40-50	20-30	30-40	40-50
















				◆											
◆		◆		◆	◆										
◆	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆◆	◆◆		◆◆		
◆	◆	◆	◆		◆	◆	◆	◆	◆	◆◆	◆◆	◆◆	◆◆	◆◆	◆◆
◆	◆	◆	◆			◆	◆	◆	◆	◆◆	◆◆	◆◆	◆◆	◆◆	◆◆
◆	◆	◆	◆							◆◆	◆◆	◆◆	◆◆	◆◆	◆◆
							◆	◆	◆	◆◆	◆◆	◆◆	◆◆	◆◆	◆◆
							◆	◆	◆	◆◆	◆◆	◆◆	◆◆	◆◆	◆◆

Lymphedema wear

“Decongestive Lymphatic Therapy (DLT) which consists of compression, exercise, manual lymphatic therapy and skin care, is the universally accepted basis of conservative treatment for lymphedema. Among its components, compression is the most important single measure which cannot be replaced by any other modality.” (Rockson, 2018).*

◆ Day Wear ◆ Night Wear

					Lymphedema		
					Mild	Moderate	Severe
					Swelling, tightness, leg heaviness, tingling in the affected extremity	Visible/Palpable skin changes, may exhibit "pitting" when pressed gently	Hard, scaly, and enlarged extremities, lymphorrhea (weeping) may be present on skin with risk of skin denudement
Product	Description		Range of Strength (mmHg)				
		Secure	Graduated compression and high containment hybrid garment	20–30	◆	◆	
				30–40	◆	◆	
				40–50	◆	◆	
		Coolflex	Graduated compression with visual indicator of compression level; breathable mesh for cool and comfortable fit	20–30		◆ ◆	◆ ◆
				30–40		◆ ◆	◆ ◆
				40–50		◆ ◆	◆ ◆
		Compreflex	Graduated, short-stretch compression; each garment includes Accutab® compression measuring devices	20–30	◆ ◆	◆ ◆	
				30–40	◆ ◆	◆ ◆	
				40–50	◆ ◆	◆ ◆	
		Comprefit	Active, inelastic, graduated compression garment, providing high levels of containment	20–50	◆ ◆	◆ ◆	
		Compreflex Reduce	Short-stretch compression garment used as a replacement for bandage wrapping; garment is customizable and used throughout therapy	Variable		◆ ◆	◆ ◆
						◆ ◆	◆ ◆
						◆ ◆	◆ ◆
		Chipsleeve®	Open-cell foam chip garment, promotes fluid movement and tissue softening, comfortable nightwear	20–25 (with Oversleeve)	◆	◆	◆
		Medaform®	Highly adjustable garment; promotes lymph drainage and provides maximum containment	30–40		◆ ◆	◆ ◆

* Stanley G. Rockson, B.-B. L. (2018). Lymphedema: A Concise Compendium of Theory and Practice (2nd ed.). [S. G. Byung-Boong Lee, Ed.] Switzerland: Springer International Publishing AG.

Custom products available upon request

Sigvaris, SIGVARIS GROUP, Accutab, Comprefit, Chipsleeve, Medaform and Compreflex are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. Coolflex is a trademark of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2023 SIGVARIS, Inc.

Rev. 06.23 LIT-LOWEXT