

# The Choice Is Clear - Veinlite<sup>®</sup> Is The Answer

Transillumination vs. Near Infrared		
	VEINLITE <sup>®</sup>	NIR PRODUCTS
<b>TECHNOLOGY USED</b>	Transillumination	NIR
<b>FIRST STICK SUCCESS RATE</b>	93%	45%
<b>WILL FIND VEINS</b>	✓	✓
<b>CLINICALLY PROVEN</b>	✓	
<b>SHOWS VEIN IN REAL TIME</b>	✓	
<b>SHOWS DEPTH OF VEIN</b>	✓	
<b>ANCHORS VEIN</b>	✓	
<b>ADDS TRACTION TO THE SKIN</b>	✓	
<b>RAISES VEIN FOR EASY ACCESS</b>	✓	
<b>PORTABLE</b>	✓	✓
<b>HAND-HELD</b>	✓	✓
<b>ONE-PERSON OPERATION</b>	✓	
<b>WARRANTY</b>	5 YEARS	1 YEAR
<b>PRICE</b>	\$499	\$3,500

Contact us to learn more about Veinlite<sup>®</sup>

281-240-3111 | [www.veinlite.com](http://www.veinlite.com) | [info@veinlite.com](mailto:info@veinlite.com)



References: A Randomized Controlled Trial. Gümüş M, Başbakkal Z. <https://www.ncbi.nlm.nih.gov/pubmed/29794951>

Peripheral intravenous cannulation with support of infrared laser vein viewing system. Rothbart A, Yu P, Müller-Lobeck L, Spies CD, Wernecke KD, Nachtigall I. <https://bmcsresnotes.biomedcentral.com/articles/10.1186/s13104-015-1431-2>