

# LabElite<sup>®</sup> DeFroster

NEW!

## Automated Rack Frost Remover

The new LabElite DeFroster removes frost from a rack that has been retrieved from long-term cold storage in either manual or automated freezers. Using mechanical brushes and a liquid solvent, the DeFroster is able to quickly and safely remove frost that may obscure barcodes on racks.

The DeFroster boasts a small footprint and is completely mobile, allowing users to perform work while using battery power. Additionally, the device's revolutionary hands-free design eliminates both the handling of harsh liquids and inhalation of noxious fumes during the defrosting process.



### Technical Specifications

Dimensions (l x w x h)	730 mm (28.7 in) x 304.9 mm (12.0 in) x 353.4 mm (13.9 in)
Supported Labware	Any SBS-format rack up to 90 mm in height
PC	Raspberry Pi
Operating System	Raspbian

### Features of the DeFroster:

- Operational using either an electrical outlet or battery power
- Won't warm samples more than 15°C, ensuring sample integrity
- Offers hands-free operation and a run time of less than 30 seconds
- The platform can operate in a cold room environment ranging from +10°C to +30°C
- Accepts SBS-format tube racks ranging in temperature from -20°C to -80°C
- Compatible with tube types up to 90 mm in height

**HAMILTON**  
STORAGE



Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)  
USA: **800-310-5866**  
Europe: **+41-58-610-10-30**

**Americas & Pacific Rim**

Hamilton Storage, Inc.  
3 Forge Parkway  
Franklin, MA 02038 USA  
Tel: +1-508-544-7000  
Fax: +1-508-544-7001  
[sales@hamilton-storage.com](mailto:sales@hamilton-storage.com)

**Europe, Middle East, Africa & Asia**

Hamilton Storage GmbH  
Parc Industrial Vial 10  
7013 Domat/Ems, Switzerland  
Tel: +41-58-610-10-30  
Fax: +41-58-610-00-30  
[sales@hamilton-storage.com](mailto:sales@hamilton-storage.com)

To find a representative in your area, please visit [hamiltoncompany.com/contacts](http://hamiltoncompany.com/contacts).