Now Buy Jordco Products Direct.



Jordco's ediore where shipping is tree.

*Free ground shipping to addresses in the USA only.

Now you can shop 24// with convenient delivery to your office. including instruments and retills, and shipping is tree. Jordco's web site. The eStore offers the tull range of Jordco products, Now you can shop for endodontic and dental products directly from Jordco is proud to announce the opening of it's new eStore.

and now online at WWW. Oraco.com Jordco products are available through your dental supply dealer




REF: ERFBHD-s (48 high density blue inserts) Jordco e-Foam Rotary HD, endodontic toam inserts, Precisely engineered to fit Jordco's Endoring organizers. tearing or shearing of the foam while a file is in motion. in cleaning debris from rotary flutes. Our Rotary HD toam minimizes paud biece can be placed directly into the high density toam to assist cleaning most rotary shaping instruments. A rotating file secured in a Jordco's e-Foam® Rotary HD inserts are specially designed for

E-FOAM ROTARY HD FREE SAMPLE

Discover for yourself how well our new e-Foam Rotary HD holds up to rotary tiles.

Email		
(^{4iD}	State	ZIP
Adress		
Name		

Fax this coupon to 503-617-1823 or scan and email to samples@jordco.com Limit one per address. Please allow up to 4-6 weeks for delivery. ~1607





To learn more about our product line, please visit www.jordco.com



ADDRESS SERVICE REQUESTED

Jordco, Incorporated

595 N.W. 167th Avenue

Beaverton OR 97006



Designed by dentists for dentists, Jordco's quality products reliably meet the challenges of your clinical practice. The Jordco System incorporates our creative and time-saving products into your current procedures. From diagnosis to treatment with safety and convenience, the Jordco System provides the ultimate organizing and storage options. Jordco's innovative dental and endodontic products are Quality Made in the USA.



ENDOGEL® FILE LUBRICATING GEL ENDOGEL lubricant is the perfect companion for the ENDORING II oraanizer. Each MicroDose tube contains 1ml of EDTA lubricating gel. ENDOGEL can be easily delivered to the pulp chamber or placed into the Gelwell cups. REF: EG-s, EG-24



TECH TIP: The Endoring II organizer vailable in Jordco blue, green, yellow and violet. Color code your code your set-ups.



ENDORING® II The ENDORING II organizer enables the clinician to directly place, store, measure and clean endodontic hand and rotary instruments within the operatina field. REF: ERK2-s (Premium Kit), ER2SK-s (Starter Kit), ER2-s (sinale with metal ruler). ER2BX-s (sinale without metal ruler)

Jordco's mission is to develop innovative dental and endodontic instruments that touch many facets of a modern practice. The Jordco System addresses most aspects of dental care and delivery. Consider adding the following Jordco products to your practice:

The Jordco System: Innovative products that will change the way you practice.

SEVENTH EDITION

ENDORING[®] MARKING TABS[™]

Three color EMTs allow you mark and identify e-Foam insert set-ups. Six numbered tabs snap off to indicate instrument usage. Each package includes 24 EMTs (8 red, 8 green and 8 blue). REF: EMT-s





TECH TIP:

The Endoring FileCaddy's 4 bulk-storage compartments allow you to organize files by size, length and type significantly reducing chair time.

TECH TIP:

Endoring Marking Tabs can be used with both the Endoring II and Endoring FileCaddy.



ENDORING® FILECADDY®

The ENDORING FILECADDY organizer provides a unique bulk storage solution for all endodontic files. Just place multiple files securely into the special foam insert, close the lid and steam autoclave the entire assembly. You're now ready for your next appointment. REF: EFC-s

To order any Jordco products, please contact your dental supply dealer. Or visit us @ www.jordco.com or call 800-752-2812

For a complete list of Jordco products with product codes, please visit www.jordco.com



own without metal rule

CONVENIENT SINGLE-USE DESIGN

control for mixing most all cements and

REF: JMS-s

JORDCO MIXING STICKS provide finger-tip

dental materials. Ideal flexibility to assure

complete mixing. 50 pack of assorted colors.

The ENDORING II organizer is available in three distinctive colors

(Yellow, Green, Violet), Color editions are available with or without a

ER2V-s (violet). REF (without metal ruler): ER2GX-s (areen). ER2YX-s (vellow).

that add the dimension of color codina to instrument set-ups

metal ruler. REF (with metal ruler): ER2G-s (areen), ER2Y-s (vellow),

ER2VX-s (violet). ER2CVP-s (areen, vellow, violet pack)

FILES INTO REPLACEABLE FOAM INSERT.

IT'S READY TO PROCESS.

CLOSE LID AND

© 2000-2016, Jardco, Incorporated. All rights reserved. Jardco products are covered by one or more U.S. or foreign patents or patent applications owned by Jardco, Inc. of Beaverton, OR, USA. See http://www.jardco.com/about.html for details. JORDCO, ENDORING, E-FOAN

ENDOGEL, E-DX, and FILECADDY are trademarks of Jordco, Inc., registered in the United States and/or other countries. The blue color of the form inserts and the endodontic organizer, the yellow color of the form inserts, and the colors yellow and blue yellow together on form insert packaging, are all trademarks of Jordco, Inc., registered in the United States. See the above website for further information. (1507) ENDORING® MARKING TABS™

EASILY IDENTIFY SET-UPS BEFORE APPOINTMENT Save time and increase efficiency. Resin designed to withstand limited steam autoclaving and allows marking with indelible ink pens.

Three color EMTs allow you mark and identify e-Foam insert set-ups. Six numbered tabs snap off to indicate instrument usage. Each package includes 24 EMTs (8 red. 8 areen and 8 blue), REF: EMT-s

ENDOGEL[®] FILE LUBRICATING GEL



IMPROVED MICRODOSE® TUBE DELIVERY SYSTEM

ENDOGEL lubricant is packaged ENDOGEL lubricant is the in 1ml MicroDose tubes to perfect companion for the reduce cross contamination . ENDORING II oraanizer. and simplify use. Each MicroDose tube

contains 1ml of EDTA lubricating gel. ENDOGEL can be easily delivered to the pulp chamber or placed into the Gelwell cups. REF: EG-s. EG-24

FEATURES	BENEFITS	
Single-Use MicroDose Tube	Reduces cross contamination	
Convenient Delivery	Polytube neck for precise placement	1
Water Soluble	Easily flushes out of canals	/

E-RULER[™] The ruler's markings are precisely laser etched providing enhanced contrast for accurate and easy calibration of endodontic hand and rotary files. REF: ERUL-s (with Stop Locks), ERULX-s (without Stop Locks).



FEATURES 30mm Laser-Etched Scale Dual-Handed Ruler

inger-Held Ring

Precise millimeter/numeric markings For use by either right or left handed operators Comfortably worn on the finaer

BENEFITS

PURE BOND DISPENSER®

Disposable, multi-purpose, light-protected storage wells for bonding agents and composites. Store lightsensitive bonding agents, sealants and composites in five wells protected by a light-locking lid. REF: PBD-s



FEATURES Increased Efficiency Light-Protecting Lid **Compact Design**

BENEFITS

Preload agent to save valuable chair and set-up time Prevents light contamination of bonding agents or composites Easily held by assistant

PURE BUFF® POLISHING WHEEL SYSTEM

The PURE BUFF polishing wheel is used for safe, efficient and cost effective light duty polishing of dental acrylics. This unique product incorporates a disposable foam wheel, autoclavable mandrel and O-ring that can be used with most dental lathes. REF: PBK-s Refill pack of 48 disposable wheels (PBF-s)

FEATURES	BENEFITS
Single-Use Wheels	Reduces cross contamination
Efficient Design	Saves valuable chair time with each adjustment
Perfect for Spot Acrylic Adjustments	Appliance can be polished without prior disinfecting

PURE FLOSS DISPENSER®

A unique floss dispensing system in which the floss is suspended between the storage and cutting towers. This suspension of floss allows the operator to use a no-touch technique to dispense floss and reduce cross contamination. REF: PFD-s



FEATURES

Durable Metal No-Touch Dispensing Sanitary Design

Zinc housing and stainless steel cutter cap Use only one hand to dispense floss

BENEFITS

Precise fit and finish reduces areas that can trap contaminates